

DRE Premier XP

Veterinary Anesthesia Machine

A reliable anesthesia machine with a small footprint

Features

- ✓ 0.1 to 4 oxygen flow meter with rotary float.
- ✓ Rotary float and large scale for easy readability.
- ✓ Recessed flush valve to prevent accidental activation.
- ✓ Quick-connect fresh gas outlet for easy and secure switching between breathing circuits.
- ✓ SELECTATEC™ compatibility manifold for easy and quick installation of vaporizer.
- ✓ Mobile five wheel spider base.
- ✓ 1500 ml CO2 absorber designed for easy and fast replenishing of sodalime. Condensation trap to prevent water build-up in sodalime. Easy and efficient gas sealing.
- ✓ Patient circuit with 2 and 3-liter bags.
- ✓ New momentary closed feature incorporated into the pop-off valve.
- ✓ Anesthesia vaporizer: New or Professionally Refurbished (Sevoflurane, Isoflurane, Halothane)

Optional Features

- ✓ Bain circuit adaptor: for non-rebreathing applications
- ✓ Scavenger interface: Active or Passive Scavenger Interface For fan or blowing assemblies
- ✓ For central vacuum system
- ✓ E-tank holder: Single and dual available
- ✓ Circuit Alarm
- ✓ Top shelf
- ✓ Wall-mount unit available
- ✓ Warranty: 5 Years
- ✓ Deeper Canister Base



SPECIFICATIONS



Weight:
12.27 kg | 27 lbs.



Length:
68.58 cm | 27"

Width:
66.04 cm | 26"

Height:
138.43 cm | 54.50"



Material:
Brass
Aluminum
Plastic
UHMW

439-2017-11-20