

### The right choice for veterinarians who need a reliable anesthesia machine with a small footprint

- » 0.1 to 4 oxygen flowmeter with rotary float.
- » Large, easy-to-read scale.
- » Recessed flush valve for preventing accidental activation.
- » Quick-connect fresh gas outlet lets you simply and securely switch between breathing circuits.
- » Selectatec™-compatible manifold makes installing the vaporizer quick and easy.
- » Mobile five-wheel base.
- » 1500 ml CO<sub>2</sub> absorber lets you easily replenish sodalime. The absorber also has a condensation trap that prevents water build-up in sodalime. It also efficiently seals gas.
- » Patient circuit with two-liter and three-liter bags.
- » Pop-off valve features new momentary closed feature.

#### Options

- » Anesthesia vaporizer
  - Isoflurane, Sevoflurane, Halothane
- » Bain circuit adaptor
- » E-tank holder
- » Circuit alarm
- » MRI compatibility
- » Wall-mount unit available



Order now to receive a  
**5-year warranty**