

Avante Vision T3

LED Surgical Light

Energy efficient technology,
ergonomic design

Mobility and Function

LED technology gives high performance lighting for the procedure area. The Vision T3's intricate design and robust materials also allow for a broad range of movement and easy positioning of the fixture.

Even Illumination of the Area

The unique optic of the Vision T3 ensures an optimum light field. The illuminated area meets the highest technical lighting requirements for treatment lights. The reflection-, shadow-, and glare-free lighting generates an ideal field of vision specific to close work.

Color Rendering

The Vision T3 offers the optimum in color rendering with Ra > 95 and R9 > 90. This allows for clear recognition of any issues involved within the tissue.

Illumination / Light Intensity Regulation

Easy operation: One control switches the light on/off and dims with five light levels, providing the practitioner with the quick, convenient control necessary when many different procedures are being done during the course of the day.

Glare Free

The light beam is focused and thus glarefree, whatever the angle of sight. This avoids all bothersome reflections on the surface.

FEATURES

- › Continuously dimmable from 100%- 10%
- › Maintenance-Free Green Technology
- › 55,000 lux/39.5in (1m)
- › Illuminated area 6.5in (160mm)
- › CRI > 95
- › Glare-free
- › Shadow and reflection-free



SPECIFICATIONS



Ceiling Mount Weight (net):
approx. 29 lbs

Floor Version Weight (net):
approx. 44 lbs



Luminous Power:
55,000 lux / 39.5 in (1 m)

Luminated Field:
6.5 in (16 cm)

Color Temperature:
approx. 4500k



Power Supply:
Direct power supply

Power Consumption:
approx. 30w



Luminaire Body:
Polycarbonate & Acrylonitrile
Butadiene Styrene
Copolymer (pc & abs)

Mounting Options:
Ceiling and Floor



620-2019-10-28