

High-Quality stainless table

This NSF approved DRE exam table has a full wrap marine edge top to help control spills, constructed of 18 gauge 304 stainless steel. It features a powder coated, weighted pedestal for attractive stability. (Constructed of 5" square, 1/4" thick tubing with a heavy-duty 14 gauge 304 stainless steel top support and a base made of extra heavy-duty 1/4" thick steel.) Includes heavy-duty adjustable feet to avoid floor damage while moving.



► A Quick Look

- » 48" long
- » 22" deep
- » 34" high

