

# Avante SLS 9000

## Portable Floor Light

### Efficient surgery lighting system provides extremely cool output

The SLS 9000 Portable LED light is energy saving and eco-friendly, with an ultra-long service life, up to 30,000 to 50,000 hours.

The CAD designed multi-LED-lighting-beam produces lighting depth more than 55 inches, with up to 90,000 lux of central illumination. New alloy construction and a simple suspension system is easy to control and be positioned accurately. This easily adjustable light system is suitable for pendant, wall-mounted and floor-standing lamps, meeting requirements of multiple operating rooms.



## SPECIFICATIONS



**Diameter:**  
20"

**Color Temperature:**  
3500°-5000° Kelvin

**LED Life:**  
50,000-hour

**Illuminance:**  
90,000 lux  
(9,000 foot candles)  
per light head

**Color Rendering Index:**  
93%

**Diameter of Light Spot:**  
5" - 14"

**LED Power:**  
55W

**Illumination Depth:**  
28" - 59"

**LED Service Life:**  
50,000 hours



*Easy-to-use  
control panel*

2196-2019-08-26