

Aeron

Veterinary Surgery Table

Equipment for the way *you* operate.

Features 20° lateral tilt and adjustable heights of up to 43"

and in the U.S. it's only available from DRE

The Pannomed Aeron Surgery Table is designed specifically for performance in the veterinary operating room environment. It has fully electric controls that allow for height adjustment up to 43". This is great for surgeons who would like to stand upright during procedures.

Additional features:

- Trendelenburg and Reverse Trend to 15° and Laterally to 20°
- X-Ray translucent top that measures 55" x 19.5" (140 cm x 50 cm)
- Operates on battery or AC power
- Braking castors provide mobility
- Angled table plane supports animals that weigh up to 250 pounds/120 kg and provides improved drainage during surgery



◀ Available with adjustable side supports



Popular with veterinary/research surgeons specializing in fractures, cardiovascular, cardiac catheterization and angiograms.

► A Quick Look

- C-Arm compatible
- Fluoroscope compatible
- 250-pound capacity
- Translucent table top
- V-top table top
- Castors / wheels
- Battery or electric options
- Lateral positioning arms