



### Reliable, MRI-compatible system is the right anesthesia choice for veterinary imaging facilities

- Constructed with non-ferrous materials
- 0.1 to 4 oxygen flowmeter with rotary float
- Large, easy-to-read scale
- Recessed flush valve for preventing accidental activation
- Quick-connect fresh gas outlet lets you simply and securely switch between breathing circuits
- Selectatec™-compatible manifold makes installing the vaporizer quick and easy
- Mobile five-wheel base
- 1500 ml CO<sub>2</sub> absorber lets you easily replenish sodalime, has a condensation trap that prevents water build-up in sodalime, and efficiently seals gas
- Patient circuit with two-liter and three-liter bags
- Pop-off valve features new momentary closed feature
- Available with MRI-compatible Penlon vaporizer
- Includes 5-year manufacturers warranty

#### Available with the Penlon Nuffield Anesthesia Ventilator

- MRI-compatible, pneumatically-driven, time-cycled ventilator
- Pre-set volume and flow rate for small or large animals
- Powered by medical air or oxygen at 340 kPa – 410 kPa



*\*Warranty does not include Nuffield Anesthesia Ventilator*