

Pannomed CCLT

Critical Care Lift-Transport Table

Equipment for the way *you* operate.

Swift and sturdy with endless accessories

The Pannomed Critical Care Lift-Transport (CCLT) is designed specifically for performance in high-traffic emergency veterinary practices. The CCLT has full electric controls that allow for height adjustment from 12" up to 39". This is great for veterinarians who need to securely transfer large animals from the floor to counter-height. Aside from flexible height adjustment, the CCLT offers castors for transport, a tilting tray for a monitor and an integrated IV pole. The synthetic plate top measures 51" by 23".



► A Quick Look

- Battery operated
- 250-lb. weight capacity
- Height adjustable from 12" to 39"
- Synthetic plate: 51" by 23"
- Lockable 5th wheel
- Tilting tray for monitor
- Integrated IV pole
- Patient loading ramp
- Patient restraint belts and tarp

Exclusively available in the U.S. from DRE!

Available with a variety of options and accessories:

- Stainless Steel Dish with Grid
- Adjustable Side Supports
- Limb Holder
- Traction Device
- Elbow Pads